




## **Technical Data**

Self-laminating labels, manual marking RiteOn, Material 1401, dispenser 22.01.2024

RO201-1401-WH

Article Number 550-14011



**Test Method** 

Material: Type 1401, Vinyl (1401)
Colour: White (WH), Transparent (CL)

Operating Temperature :  $-40 \,^{\circ}\text{C}$  to  $+80 \,^{\circ}\text{C}$  Curing Temperature: from  $+10 \,^{\circ}\text{C}$ 

Bundle ∅ min.:

Bundle Diameter max (metric):

8.10 mm

Self adhesive (Yes/No):

Yes

 Width (W):
 12.70 mm

 Height (H):
 12.70 mm

 Height (H2):
 38.10 mm

 Width of Liner (WL):
 12.70 mm

 Labels per Row:
 1 pc.

Thickness of Foil: 80 µm ASTM D3652-83

Thickness of Adhesive (TA) nom.: 27.5 µm

Adhesive: Acrylic

Initial tack: 110 g/cm<sup>2</sup> ASTM D2979-71

Peel strength: 550 N/m PSTC 7
Shear strength: 120 h FINAT FTM 7

We assume no responsibility that the Product will be used for the intended purpose, this is beyond our control. It's up to the customer to check the suitability for a specific application by means of tests. Technical changes and errors reserved





## **Technical Data**

22.01.2024

Self-laminating labels, manual marking RiteOn, Material 1401, dispenser

Print Method: manual marking

Shelf Life: Minimum 1 year upon receipt / maximal 2

years after production date.

+21 °C Storage Temperature:

Storage conditions: 50% relative humidity., The storage in the

original packaging is recommended., Please avoid warehousing under impacts such as high humidity, heat and coldness.

The functionality and durability of the labels Application Method:

can be negatively affected if improperly processed or applied. All surfaces to be bonded must be clean, dry as well as free from dust and grease. Avoid touching the adhesive surface of the label as this could

impair the application performance.

RoHS

We assume no responsibility that the Product will be used for the intended purpose, this is beyond our control. It's up to the customer to check the suitability for a specific application by means of tests. Technical changes and errors reserved